

Medical Gas Pipeline System PPM Inspection Report for:

# Test\_Facility\_Robot\_UK

22-Aug-24



# Executive Summary

BeaconMedaes Service Technician Robot\_UK Tester tested the Medical Gas Pipeline System indicated by this report in accordance with HTM02-01 Part B section 10 (Maintenance) and/or accordance with the facilities own specification (where requested) and other applicable standards or agencies. The tests performed on this system were found to be within test specifications unless otherwise noted. The results, dates and locations of the tests performed and specific data on each test can be found on the Test Summary pages.

The testing of this system may include all or part of the tests offered by BeaconMedaes. Only the tests indicated by the summary and data pages were performed on the dates indicated. The results of the tests reported represent BeaconMedaes findings on the dates and locations indicated by the report forms only. System components include:

1. Outlet : 2 pcs



ROBOT\_TAB

|                   |                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Are you a Robot ? | Added extra text on the Job dated 22-Aug-2024 101708<br><br>Dated on:22-Aug-2024 101332 -- Robot Long text testing ---- Robot Long text testing ---- Robot Long text testing ---- Robot Long text testing ---- Robot Long text testing ---<br>- Robot Long text testing ---- Robot Long text testing ---- Robot Long text testing ---- Robot Long text testing ---- Robot Long text testing -- |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|



# Inspection Report

## OUTLET

| # | Test Date | Location             | Product                                                 | Outlet condition is acceptable | Mechanical function is acceptable | Gas specificity is correct | Outlet does not leak | Flow and pressure drop is acceptable | AGSS Outlet . Flow at 1kPa | AGSS Outlet. Flow at 2kPa | AGSS Outlet. Flow at 4kPa | Risk assessment           |
|---|-----------|----------------------|---------------------------------------------------------|--------------------------------|-----------------------------------|----------------------------|----------------------|--------------------------------------|----------------------------|---------------------------|---------------------------|---------------------------|
| 1 | 22-Aug-24 | Building 1 : Floor 1 | Outlet: GEM10: Oxygen: BS: BeaconMedaes: Gem 10 T/U, O2 | Pass                           | Pass                              | Pass                       | Pass                 | Pass                                 |                            |                           |                           | I am OK to commence works |
| 2 | 22-Aug-24 | Building 1 : Floor 1 | Outlet: GEM10: Oxygen: BS: BeaconMedaes: Gem 10 T/U, O2 | Pass                           | Pass                              | Pass                       | Pass                 | Pass                                 |                            |                           |                           | I am OK to commence works |

## Description of Work

| Label                                | Description                                                                                                                                                                                                                                                                                                                                  |
|--------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Outlet condition is acceptable       | Examine terminal unit for obvious leaks or damage                                                                                                                                                                                                                                                                                            |
| Mechanical function is acceptable    | Probe operates and engages/disengages correctly and does not swivel (wall type). Outlet does not stick or jam and operates smoothly.                                                                                                                                                                                                         |
| Gas specificity is correct           | Insert a blank probe and check gas specificity is correct and no cross connections are present                                                                                                                                                                                                                                               |
| Outlet does not leak                 | Check the outlet does not leak                                                                                                                                                                                                                                                                                                               |
| Flow and pressure drop is acceptable | Using a calibrated test gun check the terminal unit performance                                                                                                                                                                                                                                                                              |
| AGSS Outlet . Flow at 1kPa           | Applicable only for AGSS outlets only                                                                                                                                                                                                                                                                                                        |
| AGSS Outlet. Flow at 2kPa            | Applicable only for AGSS outlets only                                                                                                                                                                                                                                                                                                        |
| AGSS Outlet. Flow at 4kPa            | Applicable only for AGSS outlets only                                                                                                                                                                                                                                                                                                        |
| Risk assessment                      | <p>Consider the following list as a minimum when carry out this assesment:</p> <ul style="list-style-type: none"> <li>• General Site Safety Condition</li> <li>• Hot Works</li> <li>• Housekeeping</li> <li>• Lighting</li> <li>• Noise Levels</li> <li>• Public Safety</li> <li>• Steps and Ladders</li> <li>• Working at Height</li> </ul> |



Authorization to carry out work with/without an order number.

Notes: Customer will agree to supply BeaconMedaes with a valid order number or payment for non-account holders for the actual time spend, traveling and parts used during the breakdown intervention by next working day.

**Customer Confirmation**

|

Signed by Robot\_UK  
(Being duly authorized to sign on behalf of the customer)

**Service Engineer Confirmation**

|

Signed by Robot\_UK Tester





**Test\_Facility\_Robot\_UK**

**22-Aug-2024 005838**

